





Have you heard about the P-TECH Protein Analysis?

Cert III Laboratory Skills students have! Thanks to our industry partner BioLAB our students participated in an extra special lab session this term. They acted as "scientists for the day" and they used the Bradford protein assay to measure the concentration of total protein in a sample and recorded their findings.

Thanks to mentors Jo and Craig from BIOLAB for your ongoing support of our P-TECH Program.



Using the spectrometer, find the absorbance of the standards at 595nm.

Effective Protein Concentration mg/ml	Absorbance at 595nm
0.0	
0.5	
1.0	
1.5	

